




EZ Way, Inc.

EZ Way Hip Sling Instructions

Transferring with the Hip Sling

1. Raise the bedrails on the patient's bed.
2. The sling should be positioned under the patient by log rolling the patient and positioning the sling so their hips are in the center of the sling near the center straps, and the head is inside the sling. (Figures 1 & 2)
3. Move lift to patient making sure the hangar bar is over the patient's abdomen.
4. Note: There are three attach points on each side; head, middle, and leg. Attach the corresponding loops of the hip sling to the sling hangers on the hangar bar. On a ceiling lift, the loops of the middle strap can be placed on the rear hook of the sling hanger (Figure 3). On a mobile lift, the loops of the middle strap can be placed on the sling hanger bar. (Figure 4)
5. After ensuring that all loops are connected, press the up key to raise the lift and remove slack from the straps. Once tension is placed on the straps, pause lifting and check that the loops all six straps are seated in the hooks. (Figure 5)
6. Continue to raise the patient until they have cleared the surface. Note the knees will bend. Make sure to raise the patient's feet, if needed, to help clear the surface. Maneuver the lift to complete the transfer, raising the feet, if needed to clear the surface to which you are transferring the patient. (Figure 6)



Please contact EZ Way at 1-800-627-8940 with any questions regarding sling dimensions and sizes.